

Otosillo Professional Electrisc Hoorirrigatiesysteem

Apparaat en onderdelen:

Art.nr.	Omschrijving
3039570	OTOSCILLO OORSPOEL APPARAAT PROFESSIONEEL
3039526	OTOSCILLO OORTIPS DISPOSABLE
3039527	OTOSCILLO SPATSCHERM TBV OORTIPS TRANSPARENT
3039528	OTOSCILLO HANDVATBEVESTIGING MET AUTOMATISCHE TERUGSLAG
3039529	OTOSCILLO TEMPERATUURKAART 32-41 °C
3039550	OTOSCILLO WATERCONTAINER



Disposable nozzles

According to EU recommendation
Compared to disinfecting it is more economic, **cheaper**, hygienically safer and more environmentally friendly to use a new, clean nozzle for each patient !

Automatic back stop directly behind disposable nozzle – for preventing **cross-infections** through contamination of the water-containing device system with **rinsing water of the auditory canal** which can contain bacteria, fungi and viruses !

Always the **right temperature** of the **rinsing water** with the otosillo ® Temp.-Test-Card for a quick regulation of the water temperature from the hot water tap. **Response time 0,5 sec !** For further control of water temperature attachable in the water container .

Pulsating irrigation flush by means of piston pump. Highly effective, otologically not irritating pulsation frequency via **transmission control** .

Removable water container with lid and automatic valve for water transport .

Adjustment of irrigation flush with

infinitely adjustable feed regulator for water with scaling.
Integrated purification stage for mechanical cleaning of system.
ON/OFF pump by means of **handle switch !**

Handle attachment with automatic back stop for disposable nozzles, infinitely rotatable for exact adjustment of irrigation angle. Automatic locking mechanism for nozzles !

Permanent biologic decalcification and inhibition of germ formation according to guidelines concerning hygiene .

Transparent **guard sheet**, attachable

Robust spiral water hose transparent in back light for hygiene control. Maximum extended length about 100 cm !

Guarantee 2 years – CE conformity MPG
Made in Germany

Attach



Always ! hold at the black bend of the attachment of the handle, not the handle itself ! Attachment has to be rotatable well at the colored ring. If there is any sluggishness remove attachment from handle and lubricate attachment cone at the thickened part with Vaseline ®

Remove



Pull back release ring. Remove guard sheet **with disposable nozzle**. Take nozzle at clean end and dispose it. Guard sheet can be cleaned with hot water .



Attach disposable nozzle up to the stop. Locks automatically .



Attach guard sheet **up to the stop at the release ring**. If there is any sluggishness lubricate nozzle with Vaseline ®
Do not work without guard sheet !



Medeco

Alexander Flemingstraat 2, 3261 MA Oud-Beijerland

T +31 (0)186 63 45 00 E cs.eerstelijjn@mediq.com

Direct bestellen? E-mail dan naar: bestelling@mediqmedeco.nl