

MEDECO B.V.

Product Specifications

Brand : Klinion® Protection → SENSITIVE (LONG) SENSITIVE ULTRA SAFE

Product : Nitrile, powder-free, low dermatitis potential, non-sterile, examination gloves

Packaging & Product Reference

Product Reference	Description	Size	Color	Box Dimension (mm)*	Piece(s) per Box	Outer Carton Dimension (mm)*	Box(es) per Outer Carton
102470	Klinion Personal Protection SENSITIVE Nitrile Exam Glove	XS	Indigo	235 * 110 * 75	150	390 * 230 * 245	10
102471	Klinion Personal Protection SENSITIVE Nitrile Exam Glove	S	Indigo	235 * 110 * 75	150	390 * 230 * 245	10
102472	Klinion Personal Protection SENSITIVE Nitrile Exam Glove	M	Indigo	235 * 110 * 75	150	390 * 230 * 245	10
102473	Klinion Personal Protection SENSITIVE Nitrile Exam Glove	L	Indigo	235 * 110 * 75	150	390 * 230 * 245	10
102474	Klinion Personal Protection SENSITIVE Nitrile Exam Glove	XL	Indigo	235 * 110 * 75	150	390 * 230 * 245	10
102488	Klinion Personal Protection SENSITIVE LONG Nitrile Exam Glove	XS	Indigo	270 * 125 * 80	150	420 * 260 * 280	10
102484	Klinion Personal Protection SENSITIVE LONG Nitrile Exam Glove	S	Indigo	270 * 125 * 80	150	420 * 260 * 280	10
102485	Klinion Personal Protection SENSITIVE LONG Nitrile Exam Glove	M	Indigo	270 * 125 * 80	150	420 * 260 * 280	10
102486	Klinion Personal Protection SENSITIVE LONG Nitrile Exam Glove	L	Indigo	270 * 125 * 80	150	420 * 260 * 280	10
102487	Klinion Personal Protection SENSITIVE LONG Nitrile Exam Glove	XL	Indigo	270 * 125 * 80	150	420 * 260 * 280	10
102534	Klinion Personal Protection SENSITIVE LONG Nitrile Exam Glove	XS	Blue	270 * 125 * 80	150	420 * 260 * 280	10
102535	Klinion Personal Protection SENSITIVE LONG Nitrile Exam Glove	S	Blue	270 * 125 * 80	150	420 * 260 * 280	10
102536	Klinion Personal Protection SENSITIVE LONG Nitrile Exam Glove	M	Blue	270 * 125 * 80	150	420 * 260 * 280	10
102537	Klinion Personal Protection SENSITIVE LONG Nitrile Exam Glove	L	Blue	270 * 125 * 80	150	420 * 260 * 280	10
102538	Klinion Personal Protection SENSITIVE LONG Nitrile Exam Glove	XL	Blue	270 * 125 * 80	150	420 * 260 * 280	10
165140	Klinion Personal Protection SENSITIVE ULTRA SAFE LONG Nitrile Exam Glove	XS	Blue	270 * 125 * 80	150	420 * 260 * 280	10
165141	Klinion Personal Protection SENSITIVE ULTRA SAFE LONG Nitrile Exam Glove	S	Blue	270 * 125 * 80	150	420 * 260 * 280	10
165142	Klinion Personal Protection SENSITIVE ULTRA SAFE LONG Nitrile Exam Glove	M	Blue	270 * 125 * 80	150	420 * 260 * 280	10
165143	Klinion Personal Protection SENSITIVE ULTRA SAFE LONG Nitrile Exam Glove	L	Blue	270 * 125 * 80	150	420 * 260 * 280	10
165144	Klinion Personal Protection SENSITIVE ULTRA SAFE LONG Nitrile Exam Glove	XL	Blue	270 * 125 * 80	150	420 * 260 * 280	10

*Sizes are approximate

MEDECO B.V.

Product Specifications

Product description

Disposable, latex-free, powder-free examination gloves made of nitrile butadiene rubber with no Sulphur and/or accelerators. Klinion Protection SENSITIVE nitrile accelerator-free gloves have a low-dermatitis potential known to reduce potential risk for users sensitive to chemical additives.

Intended use

The Klinion Protection examination gloves are intended to be worn on the hand for medical purposes and to help prevent contamination between patient and examiner.

User Info available at www.medeco.org

Classification

MDD 93/42/EEC Class I

PPE Regulation (EU) 2016/425 Category III

→ EU declaration of conformity at www.medeco.org

Features

- Economical and hygienic protection for hands
- Ambidextrous, suitable for left and right hand
- Finger-textured surface
- Low dermatitis potential mitigating potential risk for sensitive users
- Powder-free

Product components

- Acrylonitrile Butadiene Rubber (NBR)
- Products listed on this specification are:
 - Latex-free
 - Phthalate-free
 - Absent of human or animal tissue
 - Absent of blood or blood derivative

MEDECO B.V.

Product Specifications

Specifications

Dimensions (mm)

Width	XS	75 ± 5
	S	85 ± 5
	M	95 ± 5
	L	106 ± 5
	XL	116 ± 5
Thickness	Cuff ref	Min. 0.05
	Palm	Min. 0.06
	Finger	Min. 0.08
Length	Standard Cuff (Ref 102470 - 102474)	Min. 240
	Long Cuff (Ref 102484-102488) Long Cuff (Ref 102534-102538) Long Cuff (Ref 165140-165144)	Min. 295

Physical Properties

Unaged	Tensile Strength	Min. 14 MPa
	Elongation at Break	Min. 500 %
	Force at Break	Min. 6.0 N
Aged	Tensile Strength	Min. 14 MPa
	Elongation at Break	Min. 400 %
	Force at Break	Min. 6.0 N
Freedom from Holes	AQL 1.5	
Powder Residue	Max 2 mg/glove	

MEDECO B.V.

Product Specifications

Quality Standards

- Manufactured in accordance with Quality Assurance ISO 13485 and ISO 9001
- In compliance with:
- EN 420:2003 + A1:2009 (fit for special purpose)
 - EN ISO 374-1:2016, EN 374-2:2015, EN 374-4:2013, EN ISO 374-5:2016
 - EN 455-1:2000, EN 455-2:2015, EN455-3:2015, EN455-4:2009
 - EN 16523-1:2015

Sterilization

N/A - Product is provided non-sterile

Storage conditions

Store in a cool and dry place. Opened boxes should be kept away from fluorescent and sunlight. Gloves are packed in dispenser which is suitable for transport. Keep the gloves in the box when not in use.

Shelf life

3 years

Pictures



Ref. product 102472

Ref. product 102487

Ref. product 102535

Ref. product 165143

Approval	Name	Date	Signature
Product Manager	Levine van Buren	20-12-2018	
RA-Manager	Francis Dekker	20-12-2018	